



# PRESSALIT®

<b>Product</b> <b>Item number</b>	Pressalit Public Classic Disinfection dispenser with sensor, 700 ml, matt black  PY40725
<b>Product description</b>	VVS: 774565631 NRF: 1438615 RSK: 8920399 EAN: 5708590392769  Disinfection dispenser with stainless steel sensor for wall mounting. One-handed disinfection pump with 700 ml tank. Designed for standard hand sanitizers. Lockable solution. Including mounting kit.  Inspection opening indicates filling level.
<b>Materials</b>	European stainless steel (AISI 316L)
<b>Warranty</b> <b>Cleaning</b>	A 2-year warranty is provided for material and manufacturing defects.  This product should be cleaned with hot water or a mild cleaning solution. Do not use chlorine based or corrosive cleaners as it may result in the surface being damaged.  Subject to changes in the product range, technical changes and printing errors.

[pressalit.com](http://pressalit.com)