



PRESSALIT®

Product	Changing table mattress 1800 mm
Item number	R8704378
	VVS: 777630768 EAN: 5708590366104 HMI: 125811
Colour	Navy blue
Product description	<p>Mattresses for daycare changing tables are available in several colors and in lengths that are adapted to the different changing table lengths.</p> <p>Accessories for changing tables. Mattress 760 mm x 1800 mm x 30 mm, for 1800 mm changing table. Phthalate-free.</p>
Cleaning	<p>DCT mattresses</p> <p>Upholstery (Cover): Techmadin, phthalate-free Daily cleaning with soapy water and subsequent wiping over with a damp cloth, and wash regularly at 60-90°, without fabric softener. Shrinkage after washing may occur. Do not turn inside out before washing. Can be tumble-dried at max. 60°. Can be disinfected with antiseptic wipes. Alcohol % < 70. Do not use disinfectants intended for hard surfaces.</p> <p>Padding: (foam) The foam contains only DINP softeners which are approved by the EU for the use in articles for children. Do not wash the foam. Do not use the mattress without the cover for reasons of hygiene.</p>





PRESSALIT®

Product	Changing table mattress 1800 mm
Item number	R8704378

pressalit.com

Pressalit A/S
Pressalitvej 1
DK-8680 Ry
Tel.: +45 8788 8788
Fax: +45 8788 8789
E-mail: pressalit@pressalit.com

PRESSALIT®