




PRESSALIT®

Product Item number	Pressalit Public Premium Hygienic waste bin 4 l, brushed steel PX30925
Product description	<p>VVS:777743415 NRF:1438634 RSK:8920372 EAN:5708590392196</p>  <p>Wall-mounted sanitary bin with bag dispenser integrated on the side, removable 4 liter sanitary bin and hinged lid. Includes mounting kit and a box with plastic sanitary bags.</p>
Materials	European stainless steel (AISI 304)
Warranty Cleaning	<p>A 2-year warranty is provided for material and manufacturing defects.</p> <p>The product is to be cleaned with hot water or a mild soapy cleaning agent. Do not use chlorine based or corrosive cleaners as it may result in the surface being damaged.</p> <p>Subject to changes in the product range, technical changes and printing errors.</p>

pressalit.com

PRESSALIT®