




PRESSALIT®

Product	Wall hung toilet 540 mm
Item number	R2091000 VVS:613230000 NRF:6190522 EAN:5708590213897 HMI: 35273 NBS:N13:10,11,300
Colour	White
Product description	Toilet, wallhung, distance from wall to the front: 540 mm Excl. fixing kit
Adjustment height	
Product capacity	Max load: 400 kg.
Materials	Ceramic
Ergonomy	Wallhung toilet makes mounting possible at individual height. Leaves space for foot support, the legs of a mobile hoist, and the carer's feet.
Warranty	Pressalit makes a commitment to remedy material and manufacturing defects for a period of five years from the invoicing date.



www.pressalit.com

Pressalit A/S
Pressalitvej 1
DK-8680 Ry
Tel.: +45 8788 8788
Fax: +45 8788 8789
E-mail: pressalit@pressalit.com

PRESSALIT®