



PRESSALIT®

Product Item number	Changing table mattress 1350 mm RB1905000
Colour	White
Product description	Accessories for changing tables Mattress 760 mm x 1350 mm x 30 mm, for 1800 mm changing table with sanitary appliances Phthalate-free Accessories for changing tables Mattress 760 mm x 1350 mm x 30 mm, for 1800 mm changing table with sanitary appliances Phthalate-free



PRESSALIT®

Product	Changing table mattress 1350 mm
Item number	RB1905000
Cleaning	<p>RB1901-RB1905</p> <p>Changing table mattress, upholstery (cover) – white. 1460 PU quality (polyurethane), phthalate-free. Pressalit recommends daily cleaning with soapy water and subsequent wiping over with a damp cloth, and weekly wash at 60°, without fabric softener. Do not turn inside out before washing. Can be tumble-dried. Can be disinfected with antiseptic wipes. Do not use disinfectants intended for hard surfaces.</p> <p>Nursing table mattresses, upholstery (cover) – red or blue. NX200; a vinyl quality that is Öko-Tex 100 approved, which means that it doesn't contain phthalates or heavy metals. Pressalit recommends daily cleaning with soapy water and subsequent wiping over with a damp cloth, and weekly wash at 40°, without fabric softener. Do not turn inside out before washing. Do not tumble-dry. Can NOT be disinfected.</p> <p>Padding (foam). The foam contains only DINP softeners which are approved by the EU for the use in articles for children. Do not wash the foam. Do not use the mattress without the cover for reasons of hygiene.</p>

pressalit.com

PRESSALIT®