




PRESSALIT®

Product Item number	Modular shelves for fixed mounting, 400 mm rod RT490000
Colour	VVS:777559000 NRF:7501027 RSK:7973955 EAN:5708590154206 HMI: 124130 NBS:N13:449 White
Product description	 <p>PLUS modular shelves are practical shelves attached to a pole, they are available for fitting to PLUS wall tracks or PLUS wash basin brackets, as well as for wall-mounting.</p> <p>Modular shelf system: Rod 400 mm, for fixed mounting on wall, incl. shelf 400 mm, towel rail and tray. Incl. mounting kit for brick wall.</p>
Adjustability Product capacity	Fixed mounting Maximum load: 8 kg. Load per shelf: 4 kg. Load on towel holder: 4 kg.
Materials Features Warranty	Aluminium. Stainless steel. Laminate. The shelf and towel rail are height adjustable. Pressalit makes a commitment to remedy material and manufacturing defects for a period of five years from the invoicing date.



PRESSALIT®

Product	Modular shelves for fixed mounting, 400 mm rod
Item number	RT490000

pressalit.com

Pressalit A/S
Pressalitvej 1
DK-8680 Ry
Tel.: +45 8788 8788
Fax: +45 8788 8789
E-mail: pressalit@pressalit.com

PRESSALIT®