




PRESSALIT®

Product Item number	Modular shelves for fixed mounting, 600 mm rod RT491000
Colour	VVS:777559200 NRF:7501028 RSK:7973956 EAN:5708590154213 HMI: 124133 NBS:N13:449 White
Product description	 <p>PLUS modular shelves are practical shelves attached to a pole, they are available for fitting to PLUS wall tracks or PLUS wash basin brackets, as well as for wall-mounting.</p> <p>Modular shelf system: 600 mm, for fixed mounting on wall, incl. two shelves 400 mm, towel rail and two trays. Incl. mounting kit for brick wall.</p>
Adjustability Product capacity	Fixed mounting Maximum load: 12 kg. Load per shelf: 4 kg. Load on towel holder: 4 kg.
Materials Features Warranty	Aluminium. Stainless steel. Laminate. The shelves and towel rail are height adjustable. Pressalit makes a commitment to remedy material and manufacturing defects for a period of five years from the invoicing date.



PRESSALIT®

Product	Modular shelves for fixed mounting, 600 mm rod
Item number	RT491000

pressalit.com

Pressalit A/S
Pressalitvej 1
DK-8680 Ry
Tel.: +45 8788 8788
Fax: +45 8788 8789
E-mail: pressalit@pressalit.com

PRESSALIT®