



PRESSALIT®

Product Item number	Diaper basket, 15.8" x 8.4" x 6.9" RT767
Product description	VVS:615801767 NRF:7501034 EAN:5708590171074 HMI: 124570 NBS:N13:449 Basket for diapers, A:15.8" x B:8.4" x C: 6.9". Basket for diapers, A:15.8" x B:8.4" x C: 6.9".
Materials	Core material: Steel Surface material: Rilsan
Warranty	Pressalit makes a commitment to remedy material and manufacturing defects for a period of five years from the invoicing date.



Product	Diaper basket, 15.8" x 8.4" x 6.9"
Item number	RT767
Cleaning	<p>This product should be cleaned with hot water or a mild cleaning solution. Do not use corrosive or abrasive cleaners as it may result in the surface being damaged.</p> <p>For surface disinfection, chlorine-based disinfectants may be used according to the product specification and safety precautions for these products.</p> <p>The following cleaning agents have been tested by Pressalit and may be used, as long as they are diluted according to the instructions on the container.</p> <p>DISINFECTANTS:</p> <ul style="list-style-type: none">Clinimax Difficil-SErisan OXY+Sprit (husholdningssprit)Suma Bac D10VirkonRodalon DesinfektionssprayLysoformin SpezialCaviCideSufanios PremiumTvättssprit 70%WipeClean UniversalCleaningChlorine Disinfection WipesKlorilli (Kiilto Clean) <p>NON-DISINFECTANTS:</p> <ul style="list-style-type: none">Suma Star free D1Taski 100Taski 200Milizid Citro <p>The product must be wiped down with a soft rag wrung out in clear water immediately after cleaning, as some materials may otherwise be discolored and stained, and the surface coating may deteriorate.</p>