



PRESSALIT®

Splashguard, transparent mat, for 100 mm raised toilet seats

Item number	T9515 VVS:616810515 NRF:6190362 EAN:5708590070551	
Description	Splashguard, transparent mat, for 100 mm raised toilet seats	
	Splashguard, transparent mat, for 100 mm raised toilet seats	
Materials	Polypropylene (PP)	
Product capacity		
Weight	0,5 kg	
Gross weight	0,5 kg	
Weight per crate	0 kg	
Qty per crate	0	
Qty per pallet	0	

pressalit.com

PRESSALIT®