



Toilet seat Colani with cover With stabilising buffers and universal hinge (B83).

Item number

R36000-B83999
EAN:5708590070650
NBS:N13:300


Toilet seat, Colani, with stabilising buffers, with cover,
incl. universal hinge with 75 mm bolt.

Colour

White

Description

Toilet seat, Colani, with stabilising buffers, with cover,
incl. universal hinge with 75 mm bolt.

**Materials**

Seat/cover: Thermoset plastic (urea)

Strap: Stainless steel V2A

Product capacity

This product has been tested statically up to 375 kg. As the European standard for technical aids for disabled, BS/EN 12182:2012, requires a static test with safety factor 1.5, the nominal maximum load for this product is therefore 250 kg.

Weight

3,703 kg

Gross weight

4,103 kg

Weight per crate

20,515 kg

Qty per crate

5

Crates per pallet

12

Qty per pallet

60



PRESSALIT®

Toilet seat Colani with cover With stabilising buffers and universal hinge (B83).

Cleaning

Use a mild soap solution to clean the toilet seat. Please take care that the seat and hinges are not left damp. For disinfecting the seat, e.g. ethanol 85% can be used. Make sure that the product is wiped down with a soft rag wrung out in clear water and dried off with a soft, dry cloth immediately after disinfection. *Avoid contact of the seat and hinges with abrasive, corrosive or chlorine-based cleaners, as these may cause damage or lead to flash rust. Therefore, when cleaning the bowl, make sure that the seat and cover are kept in an upright position until all cleaning agent has been flushed away.*

Warranty

Pressalit makes a commitment to remedy material and manufacturing defects for a period of five years from the invoicing date.



www.pressalit.com

Pressalit A/S
Pressalitvej 1
DK-8680 Ry
Tel.: +45 8788 8788
Fax: +45 8788 8789
E-mail: pressalit@pressalit.com

PRESSALIT®