



PRESSALIT®


Institutional hinge

Item number

D92999
VVS: 616810934
EAN: 5708590300474

Description

Institutional hinge, stainless steel, long bolt (100 mm)



Product capacity

Weight

0,298 kg

Gross weight

0,298 kg

Weight per crate

0 kg

Qty per crate

0

Qty per pallet

0

Subject to changes in the product range, technical changes and printing errors.



www.pressalit.com

Pressalit A/S
Pressalitvej 1
DK-8680 Ry
Tel.: +45 8788 8788
Fax: +45 8788 8789
E-mail: pressalit@pressalit.com

PRESSALIT®