



# PRESSALIT®


## Institutional hinge

**Item number**

D92999  
VVS: 616810934  
EAN: 5708590300474

**Description**

Institutional hinge, stainless steel, long bolt (3.9")  
Institutional hinge, stainless steel, long bolt (3.9")

**Product capacity****Weight**

0,298 kg

**Gross weight**

0,298 kg

**Weight per crate**

0 kg

**Qty per crate**

0

**Qty per pallet**

0

Subject to changes in the product range, technical changes and printing errors.

[pressalit.com](http://pressalit.com)

# PRESSALIT®