




# PRESSALIT®

<b>Product</b> <b>Item number</b>	Freestanding column 945 mm, for fixed height PLUS support arm R9430
<b>Product description</b>	<p>VVS:615851990 NRF:6190667 EAN:5708590320427 NBS:N13:436</p>  <p>Free-standing column, grey, for fixed height PLUS support arm. Incl. mounting kit for floor.</p> <p>Free-standing column, grey, for fixed height PLUS support arm. Incl. mounting kit for floor.</p>
<b>Materials</b>	<p>Basic material: Steel</p> <p>Surface material: Rilsan</p>
<b>Features</b>	<p>Used with freestanding toilets and where it is not possible to wallmount a support arm.</p> <p>For fixed height support arm.</p>
<b>Warranty</b>	<p>Pressalit makes a commitment to remedy material and manufacturing defects for a period of five years from the invoicing date.</p>

[pressalit.com](http://pressalit.com)