




PRESSALIT®

Kipp.scharnier m/lift-off

Produktnummer

DC4999
EAN:5708590334912



Produktkapazität

Gewicht 0,14 kg

Gewicht, brutto 0,14 kg

Gewicht pro kolli 0 kg

Stück per kolli 0

Stück pro Palette 0

pressalit.com

Pressalit A/S
Pressalitvej 1
DK-8680 Ry
Tel.: +45 8788 8788
Fax: +45 8788 8789
E-mail: pressalit@pressalit.com

PRESSALIT®