




PRESSALIT®

Product Item number	Floor standing toilet with S-trap, extra high R2097000
Colour	<p>VVS:605034414 NRF:6190523 EAN:5708590335346 HMI: 105699</p> <p>White</p>
Product description	 <p>Toilet, floor-standing, extra high 460 mm. With concealed S-trap and dual flushing system preset at 3/6 litres. Can be changed to 2/4 litres Excl. mounting kit</p> <p>Toilet, floor-standing, extra high 460 mm. With concealed S-trap and dual flushing system preset at 3/6 litres. Can be changed to 2/4 litres Excl. mounting kit</p>
Materials	Ceramic
Ergonomy	A high toilet can with advantage be used by persons who would otherwise need a raised toilet seat.
Features	Extra high. With concealed S-trap and dual flushing system
Warranty	Pressalit makes a commitment to remedy material and manufacturing defects for a period of five years from the invoicing date.



PRESSALIT®

Product	Floor standing toilet with S-trap, extra high
Item number	R2097000

pressalit.com

Pressalit A/S
Pressalitvej 1
DK-8680 Ry
Tel.: +45 8788 8788
Fax: +45 8788 8789
E-mail: pressalit@pressalit.com

PRESSALIT®