



# PRESSALIT®

<b>Product</b>	Toilet paper holder with plate, polished stainless steel
<b>Item number</b>	Q10210 VVS:776482140 NRF:6190932 RSK:8744045 EAN:5708590336893  Pressalit accessories are a series of straightforward and flexible bathroom products where the functionality is the focus.
<b>Adjustment height</b>	
<b>Materials</b>	Stainless steel AISI 304, mounting knot in brass
<b>Cleaning</b>	The product is to be cleaned with hot water or a mild soapy cleaning agent. Do not use chlorine based or corrosive cleaners as it may result in the surface being damaged.



[www.pressalit.com](http://www.pressalit.com)

Pressalit A/S  
Pressalitvej 1  
DK-8680 Ry  
Tel.: +45 8788 8788  
Fax: +45 8788 8789  
E-mail: [pressalit@pressalit.com](mailto:pressalit@pressalit.com)

# PRESSALIT®