




# PRESSALIT®

<b>Product</b> <b>Item number</b>	Pressalit Choice Toilet paper holder with plate, polished steel Q10210
<b>Product description</b>	<p>           VVS:776482140            NRF:6190932            RSK:8744045            EAN:5708590336893         </p>  <p>           Let the Pressalit Choice bathroom accessory series elevate the décor of your bathroom and create coherence with a timeless look. You get a complete range of accessories to suit any bathroom.         </p>
<b>Materials</b>	Stainless steel AISI 304, mounting knot in brass
<b>Warranty</b> <b>Cleaning</b>	<p>           A 2-year warranty is provided for material and manufacturing defects.         </p> <p>           The product is to be cleaned with hot water or a mild soapy cleaning agent. Do not use chlorine based or corrosive cleaners as it may result in the surface being damaged.         </p> <p>           Subject to changes in the product range, technical changes and printing errors.         </p>

pressalit.com