



# PRESSALIT®

<b>Product</b> <b>Item number</b>	Pressalit Choice Toilet paper holder with plate, polished steel Q10210
<b>Product description</b>	VVS:776482140 NRF:6190932 RSK:8744045 EAN:5708590336893  Let the Pressalit Choice bathroom accessory series elevate the décor of your bathroom and create coherence with a timeless look. You get a complete range of accessories to suit any bathroom.
<b>Materials</b>	Stainless steel AISI 304, mounting knot in brass
<b>Warranty</b> <b>Cleaning</b>	A 2-year warranty is provided for material and manufacturing defects.  This product should be cleaned with hot water or a mild cleaning solution. Do not use chlorine based or corrosive cleaners as it may result in the surface being damaged.  Subject to changes in the product range, technical changes and printing errors.

[pressalit.com](http://pressalit.com)

## PRESSALIT®