




PRESSALIT®

Product Item number	Pressalit Choice Shelf with glassplate 58,8 x 13,5 cm, polished steel Q10400
Product description	<p>VVS:772471540 NRF:6190934 RSK:8892164 EAN:5708590336916</p>  <p>Let the Pressalit Choice bathroom accessory series elevate the décor of your bathroom and create coherence with a timeless look. You get a complete range of accessories to suit any bathroom.</p>
Materials	Stainless steel AISI 304, glass, mounting knot in brass
Warranty Cleaning	<p>A 2-year warranty is provided for material and manufacturing defects.</p> <p>The product is to be cleaned with hot water or a mild soapy cleaning agent. Do not use chlorine based or corrosive cleaners as it may result in the surface being damaged.</p> <p>Subject to changes in the product range, technical changes and printing errors.</p>

pressalit.com

Pressalit A/S
Pressalitvej 1
DK-8680 Ry
Tel.: +45 8788 8788
Fax: +45 8788 8789
E-mail: pressalit@pressalit.com

PRESSALIT®