



# PRESSALIT®

<b>Product</b> <b>Item number</b>	Pressalit Choice Toilet brush for wall with glass bowl, polished steel Q10800
<b>Product description</b>	VVS:778497140 NRF:6190942 RSK:8714259 EAN:5708590336947  Let the Pressalit Choice bathroom accessory series elevate the décor of your bathroom and create coherence with a timeless look. You get a complete range of accessories to suit any bathroom.  Elegant toilet brush for wall mounting. Opens up extra space on the floor.
<b>Materials</b>	Stainless steel AISI 304, glass, mounting knot in brass
<b>Warranty</b> <b>Cleaning</b>	A 2-year warranty is provided for material and manufacturing defects.  This product should be cleaned with hot water or a mild cleaning solution. Do not use chlorine based or corrosive cleaners as it may result in the surface being damaged.  Subject to changes in the product range, technical changes and printing errors.

[pressalit.com](http://pressalit.com)