




Product Item number	Mounting kit V8655 (4 pcs), for aerated concrete V8655
Product description	<p>EAN:5708590351391</p>  <p>Mounting kit, 4 coach screws Ø10 x 70, plugs GB14, with grey screw covers, for aerated concrete.</p> <p>Used for: MATRIX wash basin units: R4820, R4822, R4825(UK), R4827(UK) R4840, R4842, R4845(UK), R4847(UK) R4850, R4852, R4855(UK), R4857(UK) R4870, R4875(UK), R4807(UK) SELECT wash basin brackets: R4950, R4951, R4952, R4953 VALUE wash basin bracket: R1210 SELECT TL1 toilet lifters: R8012, R8022, R8212, R8222 SELECT TL2 toilet lifters: R8112, R8122, R8312, R8322 SELECT TL3 toilet lifters: R9013, R9023, R9113, R9123 SELECT TL4 toilet lifters: R9014, R9024, R9114, R9124</p> <p>Mounting kit, 4 coach screws Ø10 x 70, plugs GB14, with grey screw covers, for aerated concrete</p>
Materials	Stainless steel and plastic



PRESSALIT®

Product	Mounting kit V8655 (4 pcs), for aerated concrete
Item number	V8655

pressalit.com

Pressalit A/S
Pressalitvej 1
DK-8680 Ry
Tel.: +45 8788 8788
Fax: +45 8788 8789
E-mail: pressalit@pressalit.com

PRESSALIT®