




# PRESSALIT®

<b>Product</b>	Fixing kit V8655 (4 pcs),for aerated concrete
<b>Item number</b>	V8655
<b>Product description</b>	<p>EAN:5708590351391</p>  <p>Fixing kit, 4 coach screws Ø10 x 70, plugs GB14, with grey screw covers, for aerated concrete.</p> <p>Used for: MATRIX wash basin units: R4820, R4822, R4825(UK), R4827(UK) R4840, R4842, R4845(UK), R4847(UK) R4850, R4852, R4855(UK), R4857(UK) R4870, R4875(UK) SELECT toilet lifters: R8012, R8022, R8112, R8122 R8212, R8222, R8312, R8322 SELECT wash basin brackets: R4950, R4951, R4952, R4953 VALUE wash basin bracket: R1210</p>



[www.pressalit.com](http://www.pressalit.com)

Pressalit A/S  
Pressalitvej 1  
DK-8680 Ry  
Tel.: +45 8788 8788  
Fax: +45 8788 8789  
E-mail: [pressalit@pressalit.com](mailto:pressalit@pressalit.com)

# PRESSALIT®