



Product	Mounting kit V8655 (4 pcs), for aerated concrete
Item number	V8655
Product description	<p>EAN:5708590351391</p> <p>Mounting kit, 4 lag bolts Ø10 x 70, plugs GB14, with gray covers, for aerated concrete.</p> <p>Used for: MATRIX sink brackets: R4808, R4826, R4828 VALUE sink brackets: R1210 SELECT TL1 toiletlifters: R8032, R8042, R8232, R8242 SELECT TL2 toiletlifters: R8132, R8142, R8332, R8342 SELECT TL3 toiletlifters: R9133, R9143 SELECT TL4 toiletlifters: R9134, R9144</p> <p>Mounting kit, 4 lag bolts Ø10 x 70, plugs GB14, with gray covers, for aerated concrete</p>
Materials	Stainless steel and plastics