




PRESSALIT®

Product Item number	Mounting kit V8663 (8 pcs), for aerated concrete V8663
	EAN: 5708590351315 
Product description	Mounting kit, 8 coach screws Ø10 x 70, plugs GB14, with grey screw covers, for aerated concrete. Used for: SCT 3000 shower change tables: R8435313, R8435318 Mounting kit, 8 coach screws Ø10 x 70, plugs GB14, with grey screw covers, for aerated concrete
Materials	Stainless steel and plastic



www.pressalit.com

Pressalit A/S
Pressalitvej 1
DK-8680 Ry
Tel.: +45 8788 8788
Fax: +45 8788 8789
E-mail: pressalit@pressalit.com

PRESSALIT®