



# PRESSALIT®

|                            |                                                                                               |
|----------------------------|-----------------------------------------------------------------------------------------------|
| <b>Product</b>             | Mounting kit V8666 (8 pcs), for aerated concrete                                              |
| <b>Item number</b>         | V8666                                                                                         |
| <b>Product description</b> | EAN:5708590351285                                                                             |
|                            | Mounting kit, 8 bolts M5 x 50, plugs FTP K 6, for aerated concrete.                           |
|                            | Used for:<br>PLUS grab bars: RT018, RT024, RT036, RT042<br>PLUS modular shelves: RT490, RT491 |
|                            | Mounting kit, 8 bolts M5 x 50, plugs FTP K 6, for aerated concrete                            |
| <b>Materials</b>           | Stainless steel and plastics                                                                  |

[pressalit.com](http://pressalit.com)