




# PRESSALIT®

<b>Product</b> <b>Item number</b>	Pressalit Choice Towel holder, single, 60 x 8,5 cm, black chrome Q50500
<b>Product description</b>	VVS:775711111 NRF:6190966 RSK:8754896 EAN:5708590351759  <p>Let the Pressalit Choice bathroom accessory series elevate the décor of your bathroom and create coherence with a timeless look. You get a complete range of accessories to suit any bathroom.</p>
<b>Materials</b>	Stainless steel AISI 304, mounting knot in brass
<b>Warranty</b> <b>Cleaning</b>	A 2-year warranty is provided for material and manufacturing defects. The product is to be cleaned with hot water or a mild soapy cleaning agent. Do not use chlorine based or corrosive cleaners as it may result in the surface being damaged. Subject to changes in the product range, technical changes and printing errors.

[pressalit.com](http://pressalit.com)

Pressalit A/S  
Pressalitvej 1  
DK-8680 Ry  
Tel.: +45 8788 8788  
Fax: +45 8788 8789  
E-mail: [pressalit@pressalit.com](mailto:pressalit@pressalit.com)

# PRESSALIT®