




PRESSALIT®

Product Item number	Pressalit Choice Towel holder, single, 60 x 8,5 cm, black chrome Q50500
Product description	VVS:775711111 NRF:6190966 RSK:8754896 EAN:5708590351759  <p>Let the Pressalit Choice bathroom accessory series elevate the décor of your bathroom and create coherence with a timeless look. You get a complete range of accessories to suit any bathroom.</p>
Materials	Stainless steel AISI 304, mounting knot in brass
Warranty Cleaning	A 2-year warranty is provided for material and manufacturing defects. The product is to be cleaned with hot water or a mild soapy cleaning agent. Do not use chlorine based or corrosive cleaners as it may result in the surface being damaged. Subject to changes in the product range, technical changes and printing errors.

pressalit.com

Pressalit A/S
Pressalitvej 1
DK-8680 Ry
Tel.: +45 8788 8788
Fax: +45 8788 8789
E-mail: pressalit@pressalit.com

PRESSALIT®