




# PRESSALIT®

<b>Product</b> <b>Item number</b>	Pressalit Choice Soapdispenser, 330 ml, polished steel Q10700
<b>Product description</b>	<p>VVS:774681140 NRF:6190941 RSK:8507853 EAN:5708590356211</p>  <p>Let the Pressalit Choice bathroom accessory series elevate the décor of your bathroom and create coherence with a timeless look. You get a complete range of accessories to suit any bathroom.</p> <p>Elegant soap dispenser for wall mounting. Opens up extra space on the washbasin.</p>
<b>Materials</b>	Stainless steel AISI 304, ABS
<b>Warranty</b> <b>Cleaning</b>	<p>A 2-year warranty is provided for material and manufacturing defects.</p> <p>The product is to be cleaned with hot water or a mild soapy cleaning agent. Do not use chlorine based or corrosive cleaners as it may result in the surface being damaged.</p> <p>Subject to changes in the product range, technical changes and printing errors.</p>

[pressalit.com](http://pressalit.com)