



PRESSALIT®

Product Item number	Pressalit Choice Soapdispenser, 330 ml, brushed steel Q30700
Product description	VVS:774681100 NRF:6190957 RSK:8507854 EAN:5708590356228 Let the Pressalit Choice bathroom accessory series elevate the décor of your bathroom and create coherence with a timeless look. You get a complete range of accessories to suit any bathroom. Elegant soap dispenser for wall mounting. Opens up extra space on the sink.
Materials	Stainless steel AISI 304, ABS
Warranty Cleaning	A 2-year warranty is provided for material and manufacturing defects. This product should be cleaned with hot water or a mild cleaning solution. Do not use chlorine based or corrosive cleaners as it may result in the surface being damaged. Subject to changes in the product range, technical changes and printing errors.

pressalit.com

PRESSALIT®