



PRESSALIT®

Product Item number	Changing table mattress 1800 mm R8704000
Colour Product description	VVS: 777630760 EAN: 5708590366098 White Accessories for changing tables. Mattress 760 mm x 1800 mm x 30 mm, for 1800 mm changing table with sanitary appliances. Phthalate-free
Cleaning	R8701-R8705 Upholstery (Cover): Techmadin, phthalate-free Daily cleaning with soapy water and subsequent wiping over with a damp cloth, and weekly wash at 60°, without fabric softener. Do not turn inside out before washing. Can be tumble-dried. Can be disinfected with antiseptic wipes. Alcohol % < 50. Do not use disinfectants intended for hard surfaces. Padding: (foam) The foam contains only DINP softeners which are approved by the EU for the use in articles for children. Do not wash the foam. Do not use the mattress without the cover for reasons of hygiene.



www.pressalit.com

Pressalit A/S
Pressalitvej 1
DK-8680 Ry
Tel.: +45 8788 8788
Fax: +45 8788 8789
E-mail: pressalit@pressalit.com

PRESSALIT®