




<b>Product</b> <b>Item number</b>	Changing table mattress 1800 mm R8704000
<b>Colour</b>	VVS:777630760 EAN:5708590366098 HMI: 125811  White 
<b>Product description</b>	Mattresses for daycare changing tables are available in several colors and in lengths that are adapted to the different changing table lengths.  Accessories for changing tables. Mattress 760 mm x 1800 mm x 30 mm, for 1800 mm changing table. Phthalate-free.
<b>Cleaning</b>	R8701-R8705 <b>Upholstery (Cover):</b> Techmadin, phthalate-free Daily cleaning with soapy water and subsequent wiping over with a damp cloth, and wash regularly at 60-90°, without fabric softener. Shrinkage after washing may occur. Do not turn inside out before washing. Can be tumble-dried at max. 60°. Can be disinfected with antiseptic wipes. Alcohol % < 70. Do not use disinfectants intended for hard surfaces. <b>Padding: (foam)</b> The foam contains only DINP softeners which are approved by the EU for the use in articles for children. Do not wash the foam. Do not use the mattress without the cover for reasons of hygiene.



PRESSALIT®

<b>Product</b>	Changing table mattress 1800 mm
<b>Item number</b>	R8704000

---

[pressalit.com](http://pressalit.com)

---

Pressalit A/S  
Pressalitvej 1  
DK-8680 Ry  
Tel.: +45 8788 8788  
Fax: +45 8788 8789  
E-mail: [pressalit@pressalit.com](mailto:pressalit@pressalit.com)

PRESSALIT®