



# PRESSALIT®

<b>Product</b> <b>Item number</b>	Changing table mattress 1800 mm R8704432
<b>Colour</b> <b>Product description</b>	VVS: 777630762 EAN: 5708590366159  Medium grey  Accessories for changing tables. Mattress 760 mm x 1800 mm x 30 mm, for 1800 mm changing table with sanitary appliances. Phthalate-free
<b>Cleaning</b>	R8701-R8705 Upholstery (Cover): Techmadin, phthalate-free Daily cleaning with soapy water and subsequent wiping over with a damp cloth, and weekly wash at 60°, without fabric softener. Do not turn inside out before washing. Can be tumble-dried. Can be disinfected with antiseptic wipes. Alcohol % < 50. Do not use disinfectants intended for hard surfaces. Padding: (foam) The foam contains only DINP softeners which are approved by the EU for the use in articles for children. Do not wash the foam. Do not use the mattress without the cover for reasons of hygiene.



[www.pressalit.com](http://www.pressalit.com)

Pressalit A/S  
Pressalitvej 1  
DK-8680 Ry  
Tel.: +45 8788 8788  
Fax: +45 8788 8789  
E-mail: [pressalit@pressalit.com](mailto:pressalit@pressalit.com)

# PRESSALIT®