




PRESSALIT®

| | |
|----------------------------|--|
| Product | Soapdispenser, brushed steel |
| Item number | Q30750 |
| Product description | <p>EAN:5708590360423</p>  <p>Let the Pressalit Choice bathroom accessory series elevate the décor of your bathroom and create coherence with a timeless look. You get a complete range of accessories to suit any bathroom.</p> <p>Elegant soap dispenser for wall mounting. Opens up extra space on the washbasin.</p> |
| Materials | Stainless steel AISI 304, glass, mounting knot in brass |
| Warranty | A 2-year warranty is provided for material and manufacturing defects. |
| Cleaning | <p>The product is to be cleaned with hot water or a mild soapy cleaning agent. Do not use chlorine based or corrosive cleaners as it may result in the surface being damaged.</p> <p>Subject to changes in the product range, technical changes and printing errors.</p> |

pressalit.com

Pressalit A/S
Pressalitvej 1
DK-8680 Ry
Tel.: +45 8788 8788
Fax: +45 8788 8789
E-mail: pressalit@pressalit.com

PRESSALIT®