



PRESSALIT®

Product Item number	Pressalit Choice Bathroom set, brushed steel Q30960
Product description	<p>EAN:5708590360393</p> <p>Let the Pressalit Choice bathroom accessory series elevate the décor of your bathroom and create coherence with a timeless look. You get a complete range of accessories to suit any bathroom.</p> <p>With this bathroom set with toilet paper holder and two towel hooks you can easily create coherence in the guest toilet or in the small bathroom.</p>
Materials	Stainless steel AISI 304, mounting knot in brass
Warranty Cleaning	<p>A 2-year warranty is provided for material and manufacturing defects.</p> <p>The product is to be cleaned with hot water or a mild soapy cleaning agent. Do not use chlorine based or corrosive cleaners as it may result in the surface being damaged.</p> <p>Subject to changes in the product range, technical changes and printing errors.</p>

pressalit.com

PRESSALIT®