




PRESSALIT®

Product Item number	Safety lock for CT 4000 manually foldable change table R8486
Product description	VVS: 777545983 EAN: 5708590371191  The safety lock holds the change table securely folded against the wall. Safety lock, 13", wall-mounted, for CT 4000 manually foldable change tables
Adjustability	Fixed mounting

pressalit.com

Pressalit A/S
Pressalitvej 1
DK-8680 Ry
Tel.: +45 8788 8788
Fax: +45 8788 8789
E-mail: pressalit@pressalit.com

PRESSALIT®