



<b>Product</b>	Mattress for MCT 2 and MCT 3
<b>Item number</b>	R842717378
	VVS: 777542718 EAN: 5708590373881 HMI: 128440
<b>Colour</b>	Navy blue
<b>Product description</b>	Mattress for MCT 2 and MCT 3 change trolleys, 695 x 1680 mm. Can also be used with MSCT 1  Mattress for MCT 2 and MCT 3 change trolleys, 695 x 1680 mm.
<b>Materials</b>	The mattress cover is manufactured from 100% polyester with polyurethane, approved in accordance with Oeko-Tex 100, Class I. It is fire retardant and impermeable to dust mites. The foam filling is flexible polyester-based polyurethane foam, which is also fire retardant and approved in accordance with Oeko-tex 100, Class I.
<b>Warranty</b>	Pressalit makes a commitment to remedy material and manufacturing defects for a period of three years from the invoicing date.

