




PRESSALIT®

Product	Mattress for MCT 2 and MCT 3
Item number	R842717378
	VVS: 777542718 EAN: 5708590373881 HMI: 128440
Color	Navy blue
Product description	Mattress for MCT 2 and MCT 3 change trolleys, 27.4" x 66.1". Can also be used with MSCT 1 Mattress for MCT 2 and MCT 3 change trolleys, 27.4" x 66.1".
Materials	The mattress cover is manufactured from 100% polyester with polyurethane, approved in accordance with Oeko-Tex 100, Class I. It is fire retardant and impermeable to dust mites. The foam filling is flexible polyester-based polyurethane foam, which is also fire retardant and approved in accordance with Oeko-tex 100, Class I.
Warranty	Pressalit makes a commitment to remedy material and manufacturing defects for a period of three years from the invoicing date.



pressalit.com