




PRESSALIT®

Product	Mattress for MCT 2 and MCT 3
Item number	R842721378
	VVS: 777542728 EAN: 5708590373898 HMI: 128433
Color	Navy blue
Product description	Mattress for MCT 2 and MCT 3 change trolleys, 27.4" x 81.9". Can also be used with MSCT 1. Mattress for MCT 2 and MCT 3 change trolleys, 27.4" x 81.9".
Materials	The mattress cover is manufactured from 100% polyester with polyurethane, approved in accordance with Oeko-Tex 100, Class I. It is fire retardant and impermeable to dust mites. The foam filling is flexible polyester-based polyurethane foam, which is also fire retardant and approved in accordance with Oeko-tex 100, Class I.
Warranty	Pressalit makes a commitment to remedy material and manufacturing defects for a period of three years from the invoicing date.



pressalit.com

PRESSALIT®