




PRESSALIT®

Product	Shower mattress
Item number	R844717302
	VVS: 777542812 EAN: 5708590373904 HMI: 128437
Colour	Pewter grey
Product description	 <p>Shower mattress for MSCT 1 shower change trolley, 660 x 1660 mm. Incl. 1.6 m hose.</p> <p>Shower mattress for MSCT 1 shower change trolley, 660 x 1660 mm. Incl. 1.6 m hose.</p>
Materials	The mattress cover is welded polyester fabric with PVC-coating. It is UV-resistant and 100% impermeable to water. It contains neither CFC nor latex. The foam filling is flexible polyether-based polyurethane foam, approved in accordance with Oeko-tex 100, Class I.
Warranty	Pressalit makes a commitment to remedy material and manufacturing defects for a period of three years from the invoicing date.
Declaration of conformity	MDR-ACC This product fulfils all relevant requirements in Regulation (EU) 2017/745 on medical devices in risk class I. This product is considered to be an accessory for a medical device, as it is not itself a medical device, but is intended to be used together with several particular medical devices to specifically enable the medical devices to be used in accordance with their intended purposes.

pressalit.com