




PRESSALIT®

Product	Shower mattress for MSCT 1
Item number	R844721304
	VVS: 777542821 EAN: 5708590373935 HMI: 128434
Colour	Pearl grey
Product description	 <p>Shower mattress for MSCT 1 shower change trolley, 660 x 2060 mm. Incl. 1.6 m hose.</p> <p>Shower mattress for MSCT 1 shower change trolley, 660 x 2060 mm. Incl. 1.6 m hose.</p>
Materials	The mattress cover is welded polyester fabric with PVC-coating. It is UV-resistant and 100% impermeable to water. It contains neither CFC nor latex. The foam filling is flexible polyether-based polyurethane foam, approved in accordance with Oeko-tex 100, Class I.
Warranty	Pressalit makes a commitment to remedy material and manufacturing defects for a period of three years from the invoicing date.
Declaration of conformity	MDR-ACC This product fulfils all relevant requirements in Regulation (EU) 2017/745 on medical devices in risk class I. This product is considered to be an accessory for a medical device, as it is not itself a medical device, but is intended to be used together with several particular medical devices to specifically enable the medical devices to be used in accordance with their intended purposes.

pressalit.com

Pressalit A/S
Pressalitvej 1
DK-8680 Ry
Tel.: +45 8788 8788
Fax: +45 8788 8789
E-mail: pressalit@pressalit.com

PRESSALIT®