




# PRESSALIT®

|                            |                                                                                                                                   |
|----------------------------|-----------------------------------------------------------------------------------------------------------------------------------|
| <b>Product</b>             | Headrest for MSCT 1                                                                                                               |
| <b>Item number</b>         | R844521304                                                                                                                        |
|                            | VVS: 777542891<br>EAN: 5708590374949<br>HMI: 138482                                                                               |
| <b>Colour</b>              | Pearl grey                                                                                                                        |
| <b>Product description</b> | Headrest for MSCT 1 shower change trolley, 660 x 220 x 100 mm.<br>Headrest, 660 x 220 x 100 mm, for MSCT 1 shower change trolley. |
|                            |                                                |

[pressalit.com](http://pressalit.com)

Pressalit A/S  
Pressalitvej 1  
DK-8680 Ry  
Tel.: +45 8788 8788  
Fax: +45 8788 8789  
E-mail: [pressalit@pressalit.com](mailto:pressalit@pressalit.com)

# PRESSALIT®