




# PRESSALIT®

<b>Product</b> <b>Item number</b>	TMA 2 arm support for mounting on toilet, wider fit R1190112000
<b>Colour</b>	VVS:615724200 NRF:6190324 EAN:5708590376233 HMI: 130724  White
<b>Product description</b>	 TMA arm supports comes in various widths, and models with or without toilet seat.  TMA 2 Arm supports for mounting on toilet, incl. R37 Dania toilet seat. With BL1999 hinge. Dark grey handles.
<b>Adjustability</b>	Fixed mounting
<b>Materials</b>	Powder-coated stainless steel. Soft polyurethane. <i>Flammability-tested (ISO 8191-1). Test report available.</i> Toilet seat and cover in thermoset plastic (urea).
<b>Ergonomy</b>	The folding arms give a good support and the moulded handles are comfortable for the hands and arms.
<b>Features</b>	The support arms can be folded individually.
<b>Warranty</b>	Pressalit makes a commitment to remedy material and manufacturing defects for a period of five years from the invoicing date.
<b>Cleaning</b>	General. The product is to be cleaned with hot water or a mild soapy cleaning agent. Do not use chlorine based, corrosive or abrasive cleaners as it may result in the surface being damaged.  The product must be wiped down with a soft rag wrung out in clear water immediately after cleaning, as some materials may otherwise be discolored and stained, and the surface coating may deteriorate.
<b>Declaration of conformity</b>	CE UK MDR This product is in compliance with all relevant requirements of: - Regulation (EU) 2017/745 on medical devices. <b>For the UK:</b> - Medical Device Regulation (MDR) 2002 on medical devices.



# PRESSALIT®

<b>Product</b>	TMA 2 arm support for mounting on toilet, wider fit
<b>Item number</b>	R1190112000

---

[pressalit.com](http://pressalit.com)

---

Pressalit A/S  
Pressalitvej 1  
DK-8680 Ry  
Tel.: +45 8788 8788  
Fax: +45 8788 8789  
E-mail: [pressalit@pressalit.com](mailto:pressalit@pressalit.com)

# PRESSALIT®