




PRESSALIT®

Product Item number	TMA 2 arm support for mounting on toilet, wider fit R1191112000
Colour	VVS:615724100 NRF:6190325 EAN:5708590376240 HMI: 130727 White
Product description	 <p>TMA arm supports comes in various widths, and models with or without toilet seat.</p> <p>TMA 2 arm supports for mounting on toilet, incl. BL1999 hinge. Dark grey handles. BL1 hinge, is only recommended with toilet seats R26, R27, R28, R29, R30, R36, R37 or R38.</p>
Adjustability	Fixed mounting
Materials	Powder-coated stainless steel. Soft polyurethane. <i>Flammability-tested (ISO 8191-1). Test report available.</i>
Ergonomy	The arm supports provide firm support and the moulded handles are comfortable for the hands and arms.
Features	The arm supports can be folded up and down individually.
Warranty	Pressalit makes a commitment to remedy material and manufacturing defects for a period of five years from the invoicing date.
Cleaning	<p>General.</p> <p>The product is to be cleaned with hot water or a mild soapy cleaning agent. Do not use chlorine based, corrosive or abrasive cleaners as it may result in the surface being damaged.</p> <p>The product must be wiped down with a soft rag wrung out in clear water immediately after cleaning, as some materials may otherwise be discolored and stained, and the surface coating may deteriorate.</p>
Declaration of conformity	<p>CE UK MDR This product is in compliance with all relevant requirements of:</p> <p>- Regulation (EU) 2017/745 on medical devices.</p> <p>For the UK: - Medical Device Regulation (MDR) 2002 on medical devices.</p>



PRESSALIT®

Product	TMA 2 arm support for mounting on toilet, wider fit
Item number	R1191112000

pressalit.com

Pressalit A/S
Pressalitvej 1
DK-8680 Ry
Tel.: +45 8788 8788
Fax: +45 8788 8789
E-mail: pressalit@pressalit.com

PRESSALIT®