




<b>Product</b> <b>Item number</b>	Changing table mattress 800 mm R8700000
<b>Colour</b>	VVS:777630780 EAN:5708590376530  White 
<b>Product description</b>	Mattresses for daycare changing tables are available in several colors and in lengths that are adapted to the different changing table lengths.  Accessories for changing tables Mattress 760 mm x 800 mm x 30 mm, for 800 mm changing table Phthalate-free
<b>Cleaning</b>	R8701-R8705 <b>Upholstery (Cover):</b> Techmadin, phthalate-free Daily cleaning with soapy water and subsequent wiping over with a damp cloth, and wash regularly at 60-90°, without fabric softener. Shrinkage after washing may occur. Do not turn inside out before washing. Can be tumble-dried at max. 60°. Can be disinfected with antiseptic wipes. Alcohol % < 70. Do not use disinfectants intended for hard surfaces. <b>Padding: (foam)</b> The foam contains only DINP softeners which are approved by the EU for the use in articles for children. Do not wash the foam. Do not use the mattress without the cover for reasons of hygiene.



PRESSALIT®

<b>Product</b>	Changing table mattress 800 mm
<b>Item number</b>	R8700000

---

[pressalit.com](http://pressalit.com)

---

Pressalit A/S  
Pressalitvej 1  
DK-8680 Ry  
Tel.: +45 8788 8788  
Fax: +45 8788 8789  
E-mail: [pressalit@pressalit.com](mailto:pressalit@pressalit.com)

PRESSALIT®