




<p><b>Product</b> <b>Item number</b></p>	<p>CT 5000 change table, fixed height R8455512000</p>
<p><b>Colour</b>  <b>Product description</b></p>	<p>VVS:777540820 EAN:5708590376585 FDA:INT/D124320</p>  <p>White</p> <p>Pressalit offers a variety of mobile and wall-hung, height-adjustable and fixed height shower and changing tables of high quality. They are designed to provide a stable platform for showering and changing adults and children with reduced mobility. These tables can be used in schools, institutions, hospitals and private homes and are all designed to be durable, comfortable and easy to clean. The design and flexibility of the tables ensures optimal working conditions for the carer as well as a high level of comfort and safety for the user.</p> <p>CT 5000 change table, lying surface in compact laminate, 700 x 1200 mm , white. With safety rail and body belt, foldable, fixed height, wall-mounted.</p>
<p><b>Adjustability</b> <b>Adjustment height</b>  <b>Product capacity</b></p>	<p>Fixed mounting 0 mm</p> <p>Maximum user weight (MUM): 150 kg. Safe working load (SWL): 150 kg. When the table is loaded with 150 kg, each mounting bolt supports xx kg. <i>This product has been tested statically with a load of 225 kg as the IEC 60601-1 standard requires testing with safety factor 1.5.</i></p>
<p><b>Materials</b>  <b>Ergonomy</b>  <b>Warranty</b></p>	<p>Powder-laquered stainless steel, lying surface in compact laminate.</p> <p>The ergonomically designed safety rail has a practical lower central section, which provides better access and work environment for the carer while transferring, changing and dressing the user. Easy to clean, smooth contact surfaces which do not collect dirt and bacteria.</p> <p>Pressalit makes a commitment to remedy material and manufacturing defects for a period of three years from the invoicing date.</p>













# PRESSALIT®

<b>Product</b>	CT 5000 change table, fixed height
<b>Item number</b>	R8455512000
<b>Declaration of conformity</b>	CE UK MDR This product is in compliance with all relevant requirements of: - Regulation (EU) 2017/745 on medical devices. <b>For the UK:</b> - Medical Device Regulation (MDR) 2002 on medical devices.









<b>Product</b>	CT 5000 change table, fixed height
<b>Item number</b>	R8455512000

## Installation in different wall types

	<b>Concrete &gt;B25</b> V8597 - Mounting kit, 2 coach screws Ø10 x 70 and 2 countersunk coach screws Ø10 x 80	
	<b>Solid Brick &gt;= Mz12</b> V8597 - Mounting kit, 2 coach screws Ø10 x 70 and 2 countersunk coach screws Ø10 x 80	
	<b>Solid Sand Lime Brick &gt;= Hlz12</b> V8597 - Mounting kit, 2 coach screws Ø10 x 70 and 2 countersunk coach screws Ø10 x 80	
	<b>Aerated Concrete P4</b> V8597 - Mounting kit, 2 coach screws Ø10 x 70 and 2 countersunk coach screws Ø10 x 80	
	<b>Drywall (2 x 12,5 mm) water-resistant, re-enforced with wood</b> V8597 - Mounting kit, 2 coach screws Ø10 x 70 and 2 countersunk coach screws Ø10 x 80	



<b>Product</b> <b>Item number</b>	CT 5000 change table, fixed height R8455512000
	<b>Drywall (2 x 12,5 mm) water-resistant, re-enforced with 22 mm plywood</b> Mounting in this wall type is not recommended - risk of fatal injury to persons and walls. 
	<b>Drywall (2 x 12,5 mm) water-resistant, re-enforced with steel plates (2 mm)</b> Mounting in this wall type is not recommended - risk of fatal injury to persons and walls. 
	<b>Drywall (2 x 12,5 mm) water-resistant</b> Mounting in this wall type is not recommended - risk of fatal injury to persons and walls. 

Bolts for wall fixing are included. It is the responsibility of the installer to ensure that suitable plugs are used. Please see mounting instruction.