




PRESSALIT®

Product	Mattress
Item number	R8481378
	<p>VVS: 777540928</p> <p>EAN: 5708590376592</p> <p>HMI: 138102</p>
Color	Navy blue
Product description	 <p>Mattress for change table, 26.8" - 43.3". The mattress is easily removed, which offers flexibility of use and ensures optimum opportunity for effective cleaning over and under the mattress.</p> <p>The mattress scan be turned over and is made from a pleasant hypoallergenic material with a surface that is easy to wipe down and keep clean.</p> <p>Mattress for change table, 26.8" - 43.3". The mattress is easily removed, which offers flexibility of use and ensures optimum opportunity for effective cleaning over and under the mattress.</p> <p>The mattress scan be turned over and is made from a pleasant hypoallergenic material with a surface that is easy to wipe down and keep clean.</p>
Materials	<p>The mattress cover is manufactured from 100% polyester with polyurethane, approved in accordance with Oeko-Tex 100, Class I. It is fire retardant and impermeable to dust mites. The foam filling is flexible polyester-based polyurethane foam, which is also fire retardant and approved in accordance with Oeko-tex 100, Class I.</p>



PRESSALIT®

Product	Mattress
Item number	R8481378

pressalit.com

Pressalit A/S
Pressalitvej 1
DK-8680 Ry
Tel.: +45 8788 8788
Fax: +45 8788 8789
E-mail: pressalit@pressalit.com

PRESSALIT®