




PRESSALIT®

Product	Wall-mounted toilet, rimless, 21.3"
Item number	R2098000
	VVS:613232000 NRF:6191063 EAN:5708590377155
Color	White
Product description	 <p>Toilet, wall-mounted, distance from wall to the front: 21.3" Excl. mounting kit. For mounting directly on wall, a mounting kit from V&B must be ordered from Pressalit at an additional cost. For mounting on built-in cistern or Pressalit toilet lifter, mounting kit is included in the delivery.</p> <p>Toilet, wall-mounted, distance from wall to the front: 21.3". Excl. mounting kit.</p>
Product capacity	Max load: 880 lbs.
Materials	Ceramic
Ergonomy	Wall-mounted toilet makes mounting possible at individual height. Leaves space for foot support, the legs of a floor lift, and the caregiver's feet.
Warranty	Pressalit makes a commitment to remedy material and manufacturing defects for a period of five years from the invoicing date.

pressalit.com

PRESSALIT®