




PRESSALIT®

Product Item number	Pressalit Style Toilet paper holder, chrome PR20200
Product description	<p>VVS:776478004 NRF:7001522 RSK:8744200 EAN:5708590375137 DB: 2266842</p>  <p>Let elegance and your personal style shine through with Pressalit Style. The long, elegant, graphic shapes will lend an exciting touch to your bathroom. Not only will your bathroom get an aesthetic makeover - it will also reflect the fact that you have made an active decision to welcome extraordinary Danish design into the décor of your home. The products are easy to fit - either with the screws included or with glue, which is purchased separately.</p>
Materials	Zamak
Warranty Cleaning	<p>A 2-year warranty is provided for material and manufacturing defects.</p> <p>The product is to be cleaned with hot water or a mild soapy cleaning agent. Do not use chlorine based or corrosive cleaners as it may result in the surface being damaged.</p> <p>Subject to changes in the product range, technical changes and printing errors.</p>

pressalit.com

Pressalit A/S
Pressalitvej 1
DK-8680 Ry
Tel.: +45 8788 8788
Fax: +45 8788 8789
E-mail: pressalit@pressalit.com

PRESSALIT®