




PRESSALIT®

Product Item number	TMA 3 arm support for mounting on toilet (no call button) R1192112
Colour Product description	<p>VVS: 615725112 EAN: 5708590387475 HMI: 139267</p> <p>Anthracite grey</p>  <p>TMA arm supports comes in various widths, and models with or without toilet seat.</p> <p>TMA 3 arm supports for mounting on toilet. Anthracite grey handles. Incl. two paper holders. Excl. hinges. Bottom mounting. Compatible with IFØ Spira toilet 6261, Spira 6293 and Sign 6893. Is recommended with Pressalit Sign 553 + 554 toilet seats.</p>
Adjustability Product capacity	<p>Fixed mounting</p> <p>Maximum user weight, total: 150 kg. Maximum load per arm support: 75 kg.</p>
Materials Ergonomy Features	<p>Anodised aluminium and ABS.</p> <p>The arm supports provide firm support and the handles are comfortable for the hands and arms.</p> <p>The arm supports are foldable. The arm supports can be folded up and down individually.</p>
Warranty Cleaning	<p>Pressalit makes a commitment to remedy material and manufacturing defects for a period of five years from the invoicing date.</p> <p>General. The product is to be cleaned with hot water or a mild soapy cleaning agent. Do not use chlorine based, corrosive or abrasive cleaners as it may result in the surface being damaged.</p> <p>The product must be wiped down with a soft rag wrung out in clear water immediately after cleaning, as some materials may otherwise be discolored and stained, and the surface coating may deteriorate.</p>



PRESSALIT®

Product	TMA 3 arm support for mounting on toilet (no call button)
Item number	R1192112
Declaration of conformity	CE MDR This product is in compliance with all relevant requirements of: - Regulation (EU) 2017/745 on medical devices.

pressalit.com

Pressalit A/S
Pressalitvej 1
DK-8680 Ry
Tel.: +45 8788 8788
Fax: +45 8788 8789
E-mail: pressalit@pressalit.com

PRESSALIT®