






Pressalit 2000 without cover 49
Standard toilet seat without cover incl. hinge in stainless steel

| | |
|--|---|
| <p>Item number</p> <p>Colour</p> <p>Description</p> | <p>49000-UN3999 EAN:5708590299136 Pressalit 2000 - without cover. Light and elegant design.</p> <p>White</p> <p>PRESSALIT toilet seat without cover, model "2000", art. no. 49 made of colour ingrained duroplast with UN3 univertical hinge of stainless steel.</p> <ul style="list-style-type: none"> - centre distance: 155 mm - load - ring seat 240 kg - suitable for many universal toilets - 10 year guarantee <div style="display: flex; justify-content: space-around; align-items: flex-end;">    </div> |
| <p>Designer</p> <p>Design year</p> | <p>Pressalit</p> <p>1991</p> |
| <p>Materials</p> <p>Product capacity</p> | <p>TOILET SEAT: The material is colour ingrained duroplast (UF A 10 = urea formaldehyde) that contains no environmentally hazardous substances. UF plastic comprises 67% urea formaldehyde resin, 28% cellulose, as well as 5% minerals, pigments, lubricants and moisture content. The raw material from the supplier contains up to 0.03% free formaldehyde and even less following curing.</p> <p>BUFFERS: EVA (copolymer made of ethylene and vinyl acetate).</p> <p>DAMPER KIT (only relevant for toilet seats with soft close): Hydraulic damper with stainless steel damper casing and thermoplastic parts.</p> <p>MOUNT/HINGES: Stainless steel model W. no. 1.4301/M and thermoplastic.</p> <p>Load - ring seat 240 kg</p> |
| <p>Weight</p> | <p>1,5069 kg</p> |





PRESSALIT®

Pressalit 2000 without cover 49
Standard toilet seat without cover incl. hinge in stainless steel

| | |
|--------------------------|--|
| Carton size | 457x389x66 mm |
| Gross weight | 1,9069 kg |
| Crate size | 469x401x348 mm |
| Weight per crate | 9,5345 kg |
| Qty per crate | 5 |
| Crates per pallet | 18 |
| Qty per pallet | 90 |
| Cleaning | Use a mild soap solution to clean the seat. Please take care that the seat and hinges are not left damp. Dry off any residual water or cleaning agent with a soft cloth. Avoid contact of the seat and hinges with abrasive, corrosive or chlorine-based cleaners, as these may cause damage or lead to flash rust. Therefore, when cleaning the bowl, make sure that the seat and cover are kept in an upright position until all cleaning agent has been flushed away. |
| Warranty | <p>Pressalit provides a 10-year product warranty on your toilet seat. Pressalit will repair or exchange the toilet seat for the same or another useable model, at Pressalit's discretion, in the event of a material or manufacturing defect in your toilet seat appearing within 10 years of the date of purchase with normal use. This product warranty is conditional upon proof of purchase of the toilet seat.</p> <p>Subject to changes in the product range, technical changes and printing errors.</p> |



www.pressalit.com

Pressalit A/S
Pressalitvej 1
DK-8680 Ry
Tel.: +45 8788 8788
Fax: +45 8788 8789
E-mail: pressalit@pressalit.com

PRESSALIT®